



AVAILABLE
NOW!

gender 
SELECTed[™]
SexedULTRA[™] **4M**

New sexed-semen product available from Select Sires.

Learn more through the answers to frequently asked questions below or contact your representative today!

Is a 4 million sperm dose better than a conventional dose of Select Sires semen? No. Our study has shown that a 4 million sperm dose in Holstein heifers has a conception rate that was approximately 80 percent of the conception rate of Select Sires' conventionally processed semen. These results were statistically significant.

Does a 4 million sperm dose provide better fertility than a 2 million dose? Yes. SexedULTRA 4M has shown a conception rate advantage of 5 to 6 percentage points in Holstein heifers based on the results of a recent trial conducted by Select Sires and Sexing Technologies within some of the partner herds participating in the Program for Fertility Advancement[™].

Can this product be used on cows? SexedULTRA and SexedULTRA 4M can be used on cows. You can expect to see reduced conception rates in cows, however the results between 2 million and 4 million will be closer than in heifers. We recommend use for first or second services on cows showing strong, healthy heats.

How will this product be labeled? The '509' stud code prefix will be included for Select Sires bulls and the '559' stud code prefix for GenerVations sires. Straws and cane tops will be yellow.

Will bulls have both 4 million and 2 million available? No. At this point, bulls will either be available as SexedULTRA or SexedULTRA 4M, not both.

Now available *gender* SELECTed SexedULTRA 4M

509JE1254 JX JAMMER {4}
509JE1274 AXIS

509JE1507 SIGNAL
509JE1537 WICHITA

[™] gender SELECTed is a trademark of Select Sires Inc. gender SELECTed semen processed by Sexing Technologies for Select Sires shall be used only for the single insemination of one female bovine during natural ovulation with the intent to produce single offspring unless specifically approved on an individual customer basis by Sexing Technologies in writing. This semen is produced using the proprietary technology of XY LLC and Inguran LLC, as partially represented by U.S. Patents: 6372422, 7094527, 7208265, 7371517, 6357307, 7723116, 7758811 and 7799569. Patents pending in US and foreign markets. XY[™] is a registered trademark of XY LLC. All gender SELECTed semen is processed using SexedULTRA[™] technology. [™]SexedULTRA and SexedULTRA 4M are trademarks of Inguran LLC. Product of the USA.